



Office of Management and Budget Division of Purchasing
Non-Competitive Contract Request Form

Department	Metro Police	Department Contact	Les Skaggs
Contact Email	les.skaggs@louisvilleky.gov	Contact Phone	(502)574-7036

Contract Type: check one	New	Amendment		
		Additional Funds	Time Extension	Scope
Professional Service				
Sole Source (Goods/Services)	x			
	Start	End		
Requested Contract Dates (MM/DD/YYYY)	9/26/2023	6/30/2024		

VENDOR INFORMATION

Vendor Legal Name	Foster & Freeman USA Inc.						
DBA	Foster & Freeman USA Inc.						
Point of Contact	Amy Temenak				Email	amy.temenak@fosterfreeman.com	
Street	20145 Ashbrook Place						
Suite/Floor/Apt	#190				Phone	(888)445-5048	
City	Ashburn				State	VA	Zip Code 20147
Federal Tax ID#				SSN# (If sole proprietor)			
Louisville Revenue Commission Account #							
<u>Human Relations Commission Certified Vendors</u>				Certified Minority Owned Business	Certified Woman Owned Business	Disabled Owned Business	
Select if applicable							

FINANCIAL INFORMATION

Not to Exceed Contract Amount	45115.08		(Including reimbursement expenses, if applicable)		
Spend Category:		Law Enforcement Equipment			
Fund:	1101	Cost Center:	305 Forensics		
Program:			Project:		
Grant:			Other Worktags:		
Payment Rate		per hour		per day	X Other - describe below
		per month		per service	
Payment Frequency		Monthly	X	Upon Completion / Delivery	
		Quarterly		Other	



Office of Management and Budget Division of Purchasing

Non-Competitive Contract Request Form

CONTRACT SCOPE and PURPOSE (Attach additional documentation if necessary)

New: Be specific about the work to be performed / product to be purchased including but not limited to scope of work; description of service; work product created; why the service / product is necessary; and benefit to Louisville Metro Government.

Amendments: Describe the circumstances under which a time extension or scope change is needed.

Forensic light sources increase the efficiency and confidence of evidence collection by revealing details not visible to the naked eye. The ever-evolving field of Forensics has proven via numerous studies over the last few decades how important the use of forensic lights sources are during the search of crime scene using ultraviolet and infrared wavebands. The wavelengths can make evidence such as fingerprints, blood on dark backgrounds, other bodily fluids (semen, saliva, urine), bruising (even after the victim has healed), gunshot residue on clothing, and tattoos on decomposed decedents viewable in order to be collected and preserved for further testing and apprehension of a suspect. Foster and Freeman's portable Crime-Lite AUTO is an all-in-one imaging system that is easy to use. It is compact and easily transportable from laboratory use, to crime scenes, to autopsies. This handheld forensic light source combines the capabilities of a DSLR camera with a Crime-Lite level narrow band illumination. The Crime-Lite AUTO enables the technician to detect evidence via a live video search of the scene and capture images of evidence under multiple modes of illumination. These images are filtered to a high-quality camera and are immediately visible on the 7-inch touchscreen display with easy to use touch controls. The images can be stored on an SD card, which can then be uploaded into Axon. The use of the Crime-Lite AUTO not only accelerates the scene processing time by visually looking for bodily fluids, but also preserves the original integrity of the evidence by not contaminating it with latent powders and chemicals such as Blue Star. It can be used in light or dark environments, which is helpful in areas we cannot get pitch dark, providing a fluorescence image to enable proper evidence collection and later presentation to a jury.

JUSTIFICATION FOR NON-COMPETITIVE GOOD/SERVICE (Attach additional documentation if necessary)

Provide justification including but not limited to a description of the unique features that prohibit competition; research conducted to verify the vendor as the only known source (sole source); why the service (PSC) is not feasible to be provided by LMG staff or expertise does not exist; known compatibility, proprietary and/or timing issues.

Foster & Freeman USA Inc. is the direct manufacturer for the Crime-Lite Auto Camera.

The Crime-Lite AUTO is the most complete and compact forensics digital camera solution in existence. It combines the latest forensic imaging technology with high-intensity multi-spectral illumination, the new Crime-lite AUTO is a complete solution to the search, detection, and capture of evidence. It is a compact handheld device, with the weight and feel of a pro-grade digital SLR camera, the Crime-lite AUTO can reveal evidence including, body fluids (semen, saliva, and urine, etc.); blood on dark fabrics; fingerprints; gunshot residues; and traces of physical evidence (glass, fibers, hairs, etc.). It is capable of documenting evidence at the crime scene using a high-sensitivity UV/Vis/IR camera capable of capturing images beyond the visible spectrum. It is able to show this evidence WITHOUT using powders or chemicals to enhance them, thus it does not compromise the integrity of the evidence or the scene. The Crime-Lite AUTO uses a regular SD card that can then be used to upload photographs (or even videos) straight to AXON and then can be put under the individual case the photographs are associated with. Justification continued on next page...

AUTHORIZATIONS: Per KRS 45A.380, I have determined that competition is not feasible for the above described good / service and there is a single source within a reasonable geographical area of the good / service to be procured; or the resulting contract is for the services of a licensed professional, technician, artist, or other non-licensed professional service.

Department Director

Signature

Printed Name

Date

10/05/2023

Chief Jacquelyn Guzman-Daniel

Purchasing Director

Signature

Joel Neaveill

Date

Justification (continued)


This is the ONLY camera that allows us to use it for MULTIPLE functions. We can use it both in the lab AND on crime scenes to search, detect, AND capture evidence on site and with one device. It will allow us to only have to carry this one portable device to scenes instead of UV lights, IR lights, goggles, and our regular cameras. There is no other camera on the market that will allow us to do all of this with only one camera. We would need multiple cameras and multiple different light sources in order to do ALL of the things this one camera is capable of doing. And yes, it is the only camera that can detect saliva, blood, and even gun shot residue on dark fabrics. We can even use it to detect bruising on victims' skin that is not visible to the naked eye. It is the only camera currently on the market that allows us to do all of these things with just one camera.

foster+freeman

Crime-lite

AUTO


FOR CRIME SCENE & LABORATORY



The most integrated, complete and compact forensics
digital camera solution in existence

Crime-lite AUTO

Forensic Search and Capture



Full Spectrum UV-Vis-IR Camera with 55x Illumination Wavebands

Introducing the first in a new generation of forensic search and detection tools.

Combining the latest forensic imaging technology with high-intensity multi-spectral illumination, the new Crime-lite AUTO is a complete solution to the search, detection and capture of evidence.

A compact handheld device, with the weight and feel of a pro-grade digital SLR camera, the Crime-lite AUTO can reveal evidence including, body fluids (semen, saliva, and urine etc.); blood on dark fabrics; fingerprints; gunshot residues; and traces of physical evidence (glass, fibres, hairs, etc.)

■ Precision Imaging Camera

Document evidence at the crime scene using a UV-IR sensitive camera with autofocus macro lens.

■ UV-Vis-IR Illumination

Utilise up to 55 waveband combinations to detect a wide range of evidence types.


Search beyond the visible spectrum using UV, Vis, and IR wide-angle lighting.

■ Intelligent Interface

View crystal clear live video and captured images on the Crime-lite AUTO's simply controlled, integral 7" touchscreen interface.




Detect, Examine, and Capture Evidence at the Crime Scene



DETECT BODY FLUIDS Fluorescence Examination

Powerful multispectral illumination reveals naturally fluorescent body fluids including semen, saliva, urine, vaginal fluids etc.



REVEAL BLOOD STAINS Live Infrared Examination

Realtime infrared examination can quickly locate the presence of blood against dark backgrounds.



VISUALISE FINGERMARKS Contrast Examination

Powerful illumination and camera technology enables the examiner to visualise fingerprints on a variety of substrates.

Reveal and Enhance the Visibility of Evidence

Performing real-time visible, fluorescent, and infrared examinations, the Crime-lite AUTO can quickly reveal the presence of evidence including body fluids (semen, saliva, sweat etc.), fingerprints, fragments of bone and tooth, hair, and suspect fibres, etc.

Optimised for mobile use, to be operated by examiners in full CSI-PPE kit, the compact device further improves the evidence workflow and increases the efficiency of examinations, by allowing video and images to be stored internally or exported to an external device for immediate digital transfer to the forensic laboratory.




- ✓ Body Fluids
- ✓ Fingerprints
- ✓ Trace Evidence

Powerful Forensic DSLR Camera | Multi-Spectral Light Source

Wide-Angle Illumination
with UV, White and IR LEDs

Multi-Spectral Illumination
provides up to 55 waveband combinations

Precision 20MP Camera
sensitive from UV-IR



Crime Scene Photography Without Limitations

Search for traces of evidence under UV, Visible and Infrared illumination, before documenting the scene using multi-spectral of photography to reveal evidence that cannot be seen by the naked eye.

Optimise Illumination to Reveal More Evidence


Select LEDs individually or in combination to create up to 55x wavebands of illumination for the detection of evidence, and for the excitation of fingerprint powders and naturally fluorescent body fluids.

Integral Touch Screen Display | Advanced Connectivity

Advanced Connectivity
Bluetooth, Micro SD, USB & HDMI


Motorised Filter Wheel
with intelligent filter selection

Integrated 7" Multi-Touch Display
with forensic capture software



Intelligent Fluorescence Filter Selection

When alternating between illumination wavebands, the Crime-lite AUTO automatically selects the corresponding fluorescence imaging filter to optimise evidence visualisation. Full manual control of filter selection is also available.



Uncomplicated 7" multi-touch interface

Intuitive and easy to use, the Crime-lite AUTO's large touch screen interface provides complete control of all camera and illumination functions and settings whilst displaying crystal clear live and stored images of evidence.

Crime-lite AUTO System Specifications

System Hardware

CAMERA

- High sensitivity colour/IR camera
- Autofocus macro lens
- High resolution still images (5456x3632px)
- HDMI live video output

ILLUMINATION

High intensity LED illumination ring:

- UV 365nm, Violet 410nm, Blue 445nm, Blue Green 475nm, Green 520nm, Orange 590nm, Red 640nm
- 4x LED per wavelength
- Combine LEDs to provide up to 55 combinations
- Brightness control function

2x Wide-angle LED illumination arms:

- UV 365nm, White VIS, IR 850nm
- 2x UV LED, 1x White LED, and 1xIR LED per arm

IMAGING FILTERS

9-position motorised filter wheel

- Longpass filters at:
VIS, 455nm, 495nm, 550nm, 590nm, and 780nm
- Bandpass filters at:
UV 315-380nm, 415nm, and 530nm

INTERFACE

7" multi-touch IPS TFT display
Internal processor (no PC required)
App' style user interface
WiFi and Bluetooth connectivity

Dimensions



- H 126mm x W 207mm x D 159mm (excluding battery)
- Weight (when fitted with 2Ah battery) 1.6kg (approx.)

System Software

Crime-lite AUTO functions are carried out via an 'App' style interface using the system's integrated multi-touch display:

- Live video search and capture (AVI)
- High-resolution capture of still images (BMP, TIFF, PNG, JPG)
- Complete control of all camera settings
- Complete control of all illumination settings
- Manual or semi-automated 'intelligent' filter selection
- Pull-out image gallery of all recent work
- Save and Recall images with examination settings
- Timestamp captured video and images
- Export images and data to external USB
- Option to connect secondary viewing monitor via HDMI

Optional Accessories

SPARE 2AH BATTERY

High capacity 12V V 2Ah Lithium Ion battery with battery fuel gauge and battery protection circuit

QCL/AUTO/ACC/1

SPARE 4AH BATTERY

High capacity 12V 4Ah Lithium Ion battery with battery fuel gauge and battery protection circuit

QCL/AUTO/ACC/2

RAPID BATTERY CHARGER

Charge Crime-lite AUTO Li-ion batteries. Includes charge status display

QCL/AUTO/ACC/3

AC MAINS POWER ADAPTOR

Operate the Crime-lite auto via AC mains power

QCL/AUTO/ACC/8

CARRY CASE

Rugged, waterproof case with machined foam cut-outs

QCL/AUTO/ACC/4

ADJUSTABLE SHOULDER/NECK STRAP

QCL/AUTO/ACC/5

HIGH SPEED 128GB MICRO SD CARD

Store >2000 uncompressed images

QCL/AUTO/ACC/6

HIGH SPEED USB MEMORY STICK

QCL/AUTO/ACC/7



Head Office, UK Sales Office
Vale Park | Evesham
WR11 3TD | United Kingdom

T: +44 (0)1386 768 050
E: sales@fosterfreeman.com

USA Sales Office
40030 Manekin Plaza | Suite 170
Sterling | VA 20166 | USA

T: 888 445 5048
E: usoffice@fosterfreeman.com

EU Sales Office
Foster+Freeman GmbH
Düsseldorf | Germany

T: +49 (0)211 358 422 57
E: owen.lang@fosterfreeman.com

foster+freeman
Forensic Science Innovation

foster+freeman
FORENSIC SCIENCE INNOVATION

FOSTER+FREEMAN USA INC.
20145 ASHBROOK PLACE #180
ASHBURN, VA 20147

T: 888 445 5048
F: 888 445 5049

usoffice@fosterfreeman.com

Federal Tax ID No. 65-2318854

GST/HST No. 81959-0300 RT0001

www.fosterfreeman.com

Your Reference: RFQ: CRIME-LITE AUTO KIT

Our Reference: C103516

Quote Reference: Q214475

Date: August 29, 2023

Page 1 of 9

Louisville Metro Police Dept.
701 West Ormsby Avenue
Louisville
KY 40203
USA

Attention of Jacob Brandt

QUOTATION Q214475 - FOR THE SUPPLY OF:

LINE	PRODUCT DESCRIPTION	QTY	UNIT PRICE	SUB TOTAL
1	CRIME-LITE AUTO KIT WITH OBLIQUE AND COAX 120V (QCL/AUTO/KIT/SOC120V)	1	42,815.08	42,815.08

**CRIME-LITE AUTO KIT 120V CHARGER
(QCL/AUTO/KIT120V)**

Comprising:

**CRIME-LITE AUTO
(QCL/AUTO)**

All-in-one forensic evidence imaging system for CSI and laboratory use

CAMERA AND LENSES

HIGH SENSITIVITY COLOUR/IR CAMERA

- Precision 20.1MP image capture (5456x3632px)
- Sensitive from 350 to 1000nm
- HDMI Live video output
- Autofocus 30mm f/3.5 macro lens
- Shutter speed/exposure time up to 30 seconds

Selectable advanced camera controls

- Exposure control
- Auto and Fine Manual focus modes
- White balance

ILLUMINATION AND FILTRATION

HIGH INTENSITY LED ILLUMINATION RING

- Wavelengths: UV 365nm, Violet 410nm, Blue 445nm, Blue Green 475nm, Green 520nm, Orange 590nm, Red 640nm
- 4x LED per wavelength
- Combine LEDs to provide up to 55 wavelength combinations
- Brightness control function

2X WIDE-ANGLE LED ILLUMINATION ARMS

- Wavelengths: UV 365nm, White, IR 860nm
- 2x UV LED, 1x White LED, 1x IR LED per arm



FOSTER+FREEMAN USA INC.
20145 ASHBROOK PLACE #190
ASHBURN, VA 20147
T: 888 445 5048
F: 888 445 5049
usoffice@fosterfreeman.com

Federal Tax ID No. 56-2316854

GST/HST No. R1929 0300 RT0001

www.fosterfreeman.com

Your Reference: RFQ: CRIME-LITE AUTO KIT

Our Reference: C103516

Quote Reference: Q214475

Date: August 29, 2023

Page 2 of 9

MOTORISED FILTER WHEEL

- Motorised housing with 9-position filter wheel

FLUORESCENCE IMAGING FILTERS

- Longpass filters at: 420nm, 455nm, 495nm, 550nm, 590nm, and 780nm
- Bandpass filters at: 415nm, 530nm and UV 315-380nm

INTERFACE

TOUCHSCREEN DISPLAY

- Integrated 7" multi-touch IPS TFT display (1024x600)
- HDMI live video output
- Micro SD input, USB and Micro-USB connectivity
- Ambient light sensor with auto screen brightness
- WiFi and Bluetooth connectivity

INTERNAL PROCESSOR

- Fully-integrated processor
- Image capture (BMP, TIFF, PNG, JPG) to SD Card
- Video capture (AVI) 1920x1080 to SD Card
- Real time clock for image capture time-stamping

PHYSICAL AND ELECTRICAL

SYSTEM DIMENSIONS

- H 126mm x W 207mm x D 159mm (excluding battery)
- Weight (when fitted with 2Ah battery) 1.75kg (approx.)
- 1/4-20 UNC tripod mounting thread

ELECTRICAL

- Designed for use with 2Ah or 4Ah Li-Ion removable battery
- Up to 2-hours battery runtime
- 12Vdc power input

ACCESSORIES

- Magnetic Lens cap
- 2m Premium HDMI Lead
- UV safety glasses
- Micro fibre cloth
- Evidence scale stickers
- 128GB MICRO SD CARD SDHC/SDXC USH-I CLASS 10

SYSTEM INTERFACE


Crime-lite AUTO functions are carried out via an 'App' style interface using the system's integrated 7" multi-touch display.

OPERATING FEATURES

- Live video search and capture (AVI)
- High-resolution capture of still images (BMP, TIFF, PNG, JPG)



FOSTER+FREEMAN USA INC.
20145 ASHBROOK PLACE #100
ASHBURN, VA 20147
T: 888 445 5048
F: 888 445 5048
usoffice@fosterfreeman.com



Your Reference: RFQ: CRIME-LITE AUTO KIT

Our Reference: C103516

Quote Reference: Q214475

Date: August 29, 2023

Page 3 of 9


- Complete control of all camera settings including: Exposure, white balance, focus, shooting mode (single-shot or continuous),
- Complete control of all illumination settings including: Select single or multiple LED combinations, brightness adjust.
- Complete control of all filtration settings including: Manual or semi-automated 'intelligent' filter selection
- Pull-out image gallery of all recent work allowing comparison with current and previous images
- Save and Recall images complete with all examination settings
- Timestamp captured video and images
- Image hashing
- Option to export images and data to external USB
- Option to connect secondary viewing monitor via HDMI

PRODUCT VERIFICATION LIBRARY (PVL) 1 YEAR SUBSCRIPTION

- One-year free subscription to the foster+freeman Product Verification Library (PVL)
- Two-factor authentication for enhanced security
- Compatible with past and previous version of Chrome, Firefox, Edge, Internet Explorer, and Safari
- Over 18,000 high resolution images of body fluid stains on representative fabric types
- Over 300 high resolution sample overview (entirety) images
- 16 different garment types utilised
- 10 different garment colours of each garment type represented
- 6 different body fluids deposited (blood, saliva, semen, sweat, urine, vaginal material)
- Unique icon driven interface
- All samples and datasets compiled using Human Tissue Authority (HTA) approved body fluids
- All images captured utilising the CL-AUTO
- Multi-spectral examination of all samples recorded (UVA-NIR)
- Results analysed by forensically trained personnel to distinguish stains 'of value'
- Basic image processing utilised to represent casework workflow (greyscale and invert)
- Product Verification Library access via unique login (one user login per Crime-lite AUTO device)
- Dataset(s) can be scrutinised via simple user interface
- Results from the dataset(s) facilitate substantiated decision making at the crime scene and the laboratory
- Product Verification Library enables use of the Crime-lite AUTO to be verified by science
- Quantitative conclusions from the Ground Truth Data allow the generation of pre-sets that are scientifically valid
- Product Verification Library generated pre-sets are integrated within the Crime-lite AUTO



FOSTER+FREEMAN USA INC.
20145 ASHBROOK PLACE #180
ASHBURN, VA 20147
T: 888 445 5048
F: 888 445 5049
usoffice@fosterfreeman.com



Your Reference: RFQ: CRIME-LITE AUTO KIT

Our Reference: C103516

Quote Reference: Q214475

Date: August 29, 2023

Page 4 of 9

- Pre-sets enable the simple, efficient, and effective use of Crime-lite AUTO for the search, detection, and imaging of bodily fluids at the crime scene and in the laboratory.

LED Safety: Crime-lite AUTO is provisionally classified to European safety standard risk group 2 (IEC/EN 62471:2008)

CRIME-LITE AUTO RECHARGEABLE 2Ah BATTERY
(Quantity x 2.0)
(QCL/AUTO/ACC/1)

- 10.8 V 2Ah Lithium Ion battery
- High capacity
- Individual cell monitoring
- Battery protection circuit
- Battery fuel gauge (4-stage LEDs)

CRIME-LITE AUTO RECHARGEABLE 4Ah BATTERY
(QCL/AUTO/ACC/2)

- 10.8 V 4Ah Lithium Ion battery
- High capacity
- Individual cell monitoring
- Battery protection circuit
- Battery fuel gauge (4-stage LEDs)

CRIME-LITE AUTO CARRY CASE
(QCL/AUTO/ACC/4)

- machined foam cut-outs to accommodate Crime-lite AUTO plus batteries and accessories
- rugged, waterproof and shock resistant
- dimensions 430X260X235mm
- weight 2.7kg

CRIME-LITE AUTO SHOULDER/NECK STRAP 257-187-00
(QCL/AUTO/ACC/5)

- Adjustable shoulder/neck strap

CRIME-LITE AUTO HIGH SPEED MICRO SD CARD
READER
(QCL/AUTO/ACC/6)- Micro SD card Reader

CRIME-LITE AUTO HIGH SPEED USB STICK
(QCL/AUTO/ACC/7)- High capacity USB stick with rapid read/write performance

CRIME-LITE AUTO MAINS ADAPTOR
(QCL/AUTO/ACC/8)- Operate the Crime-Lite Auto via AC mains power

CRIME-LITE AUTO BATTERY CHARGER
(QCL/AUTO/ACC/3/120V)
- Rapid charging of Crime-lite auto Li-Ion batteries
- Charge status display



FOSTER+FREEMAN USA INC.
20145 ASHBROOK PLACE #190
ASHBURN, VA 20147
T: 888 445 5048
F: 888 445 5049
usoffice@fosterfreeman.com



Your Reference: RFQ: CRIME-LITE AUTO KIT

Our Reference: C103516

Quote Reference: Q214475

Date: August 29, 2023

Page 5 of 9

L-SHAPE FORENSIC SCALE RULE (ASSY DWG 035-106-AS-01)
(DCS/RULE)- L Shaped Forensic Rule

PELI CASE ADDRESS LABEL
(035-102-AW-01)

CRIME-LITE AUTO OBLIQUE LIGHTING RING ACCESSORY

(QCL/AUTO/ACC/8)- Oblique ring light accessory for use with CRIME-LITE AUTO

- 5700K White OR UV 365nm (Selectable – one at a time)
- 12 LEDs per wavelength
- Intensity control
- Lighting direction control – Quadrant selection
- Quadrant, wavelength and brightness controls appear on CRIME-LITE AUTO screen when connected.

CRIME-LITE AUTO COAX LIGHTING ACCESSORY

(QCL/AUTO/ACC/10)- Clip-on accessory module for Crime-lite Auto

- 5700K diffuse white coaxial light
- On-screen 'on-axis' and 'off-axis' illumination mode control (0-9)
- On-screen illumination intensity control
- On-screen focus adjustment
- Powered via USB connection to Crime-lite AUTO

INTERACTIVE EVIDENCE SAMPLE KIT

(QCL/AUTO/INT-SAMPLES) Video presentation box for Crime-lite Sample Kit.

Custom-made video box with integral 7" HD screen and two foam inserts to accommodate up to 20 evidence samples.

Video Screen:

- 7" HD display
- Pre-loaded video content
- Re-chargeable Li-Ion batteries
- Approximately 2-hours runtime

ACCELERANT (SIMULATED) (QCL/SAMPLE/A)

- Accelerant (simulated) on blue fabric to be used with blue light source (495nm LP or 530nm BP filter), blue-green light source (550nm LP or 550nm BP filter) or green light source (590nm LP filter) light source


FLUORESCENT FINGERMARK (QCL/SAMPLE/B)

- Fluorescent fingerprint (plastic sample box) to be used with UV light source (visible pass or 420nm LP filter), violet light source (455nm LP), blue light source (495nm or 530nm BP filter), blue-green light source (550nm LP or 550nm BP filter) or green light source (590nm LP filter)

POLYCYANO FINGERMARK (QCL/SAMPLE/C)



FOSTER+FREEMAN USA INC.
20145 ASHBROOK PLACE #190
ASHBURN, VA 20147
T: 888 445 5048
F: 888 445 5049
usoffice@fosterfreeman.com



Your Reference: RFQ: CRIME-LITE AUTO KIT

Our Reference: C103516

Quote Reference: Q214475

Date: August 29, 2023

Page 6 of 9

- Polycyano (plastic bag) to be used with UV light source (visible pass or 420nm LP filter)

TEXTILE FIBRES (QCL/SAMPLE/D) (QCL/SAMPLE/D)

- Textile Fibres (green fabric with fibres) to be used with violet light source (455nm LP filter), blue light source (495nm LP or 530nm BP filter) or blue-green light source (550nm LP or 550nm BP filter)

GENERAL DEBRIS (QCL/SAMPLE/E)

- Debris to be used with UV light source (visible pass or 420nm LP filter), violet light source (455nm LP filter), blue light source (495nm LP or 530nm BP filter), blue-green light source (550nm LP or 550nm BP filter) or green light source (590nm LP filter)

FOOTWEAR MARK (GEL LIFT) (QCL/SAMPLE/F)

- Footwear mark (gel lift) with white light source with visible pass filter

BODY FLUID (SIMULATED) (QCL/SAMPLE/G)

- Body fluid (simulated) on white fabric with black backing to be used with UV light source (visible pass or 420nm LP filter), violet light source (455nm LP filter), blue light source (495nm LP or 530nm BP filter) or blue-green light source (550nm LP or 550nm BP filter)

DRUGS SAMPLE (SIMULATED) (QCL/SAMPLE/H)

- Drugs (simulated) in tiny grip seal polybag to be used with UV light source (visible pass or 420nm LP filter) or violet light source (455nm LP filter)

BODY FLUID (SIMULATED) (QCL/SAMPLE/I)

- Body fluid (simulated) on black cotton to be used with UV light source (visible pass or 420nm LP filter), violet light source (455nm LP filter), blue light source (495nm LP or 530nm BP filter) or blue-green light source (550nm LP or 550nm BP filter)

BODY FLUID (SIMULATED) (QCL/SAMPLE/J)

- Body fluid (simulated) on pink and black patterned fabric to be used with UV light source (420nm LP filter), violet light source (455nm LP), blue light source (495nm or 530nm BP filter) or blue-green light source (550nm LP or 550 BP filter) – demonstrates use of bandpass for best visualisation

CYANOACRYLATE FINGERMARK (QCL/SAMPLE/K)

- Cyanoacrylate fingerprint on plastic grip-seal bag to be used with UV light source (330-385nm BP filter)

Fp NATURAL 1 FINGERMARK (QCL/SAMPLE/L)

- Fp Natural 1 fingerprint on aluminium foil to be used with blue light source (715 or 780nm LP filter) or red light source (715nm or 780nm LP filter)



FOSTER+FREEMAN USA INC.
20145 ASHBROOK PLACE #190
ASHBURN, VA 20147

T: 888 445 5048
F: 888 445 5049

usoffice@fosterfreeman.com

Federal Tax ID No. 56-2316854

GST/HST No. R1939 6300 R10001

www.fosterfreeman.com

Your Reference: RFQ: CRIME-LITE AUTO KIT

Our Reference: C103516

Quote Reference: Q214475

Date: August 29, 2023

Page 7 of 9

Fp NATURAL 2 FINGERMARK (QCL/SAMPLE/M)

- Fp Natural 2 fingerprint on wrapping paper to be used with red light source (715nm or 780nm LP filter) or IR light source (850nm filter)

NINHYDRIN FINGERMARK (SIMULATED) (QCL/SAMPLE/N)

- Ninhydrin fingerprint (simulated) on green cardboard with green light source and visible pass filter for colour-contrast imaging

BLOOD STAIN SIMULATED (QCL/SAMPLE/O)

- Blood stain (simulated) on black with red pattern fabric to be used with IR light source (IR pass filter to block out visible light either 715nm or 780nm LP filter)

GUNSHOT RESIDUE SAMPLES (SIMULATED) (QCL/SAMPLE/P)

- GSR (simulated) on black fabric with white backing to be used with IR light source (IR pass filter to block out visible light either 715nm or 780nm LP filter)

BLOOD STAIN (SIMULATED) (QCL/SAMPLE/Q)

- Blood stain (simulated) on red patterned fabric to be used with IR light source (IR pass filter to block out visible light either 715nm or 780nm LP filter)

PAINT (QCL/SAMPLE/R)

- Paint on cardboard – two paints; midnight sapphire (blue) under dark angel (purple) on top to be used with IR light source (IR pass filter to block out visible light either 715nm or 780nm LP filter)

3D FINGERMARK (OBLIQUE ACCESSORY) (QCL/SAMPLE/S)


- Oblique (3D fingerprint in mikrosil) with white light source and visible filter or UV light source (visible pass or 420nm LP filter) with oblique accessory for the CL Auto only.

RECOVER FINGERMARK (COAXIAL ACCESSORY) (QCL/SAMPLE/T)

- Coaxial (brass plate) – Recover fingerprint to be used with white light source and visible pass filter with coaxial accessory for the CL Auto only.

PRINTED INSTRUCTION MANUAL STORED IN THE SAMPLE BOX FOR EASY REFERENCE

2	PRODUCT TRAINING	1	2,050.00	2,050.00
	THIS ITEM IS EXEMPT FROM TAX			



Your Reference: RFQ: CRIME-LITE AUTO KIT

Our Reference: C103516

Quote Reference: Q214475

Date: August 29, 2023


Page 8 of 9

1-day site visit by Foster & Freeman Applications Specialist. Includes unpacking and setup of equipment, followed by comprehensive operator training and optional testing. Training will cover both theory and practical applications, using realistic samples and scenarios. Equipment safety and maintenance will also be covered. Training and testing certificates provided. Costs are based on a single visit and include associated expenses for the Applications Specialist.

3	SHIPPING AND HANDLING	1	250.00	250.00
---	-----------------------	---	--------	--------



FOSTER+FREEMAN USA INC.
20145 ASHBROOK PLACE #180
ASHBURN, VA 20147
T: 888 445 5048
F: 888 445 5049
usoffice@fosterfreeman.com



Your Reference: RFQ: CRIME-LITE AUTO KIT

Our Reference: C103516

Quote Reference: Q214475

Date: August 29, 2023

Page 9 of 9

TOTAL PRICE	USD	45,115.08
TOTAL TAX	USD	0.00
TOTAL VALUE	USD	45,115.08

DESPATCH PERIOD:

8 working weeks from receipt of order, subject to terms of payment details below. Please note despatch times are estimates and will be confirmed at time of order acknowledgement.

TERMS OF DELIVERY:

This price covers delivery including customs clearance, duties and local delivery, but excludes local sales tax.

TERMS OF PAYMENT:

Nett 30 Days from date of invoice Please note that a surcharge of 3% will apply to all credit card payments

BANK DETAILS:

United Bank, McLean, VA. Routing No 056004445, Foster & Freeman USA Inc., Account No. 70061289

WARRANTY:

24 months from date of despatch

VALIDITY UNTIL:

November 29, 2023

SIGNED:

Jason Cole

Foster & Freeman Ltd. terms and conditions of sale apply.